

A detailed view of the Hummingbird Synergy Ventricular probe. The probe is primarily white with a blue handle at the top and a blue connector at the bottom. It features a long, thin, flexible shaft with a textured section. The main body has a complex, multi-ported design with a small blue bird logo. The background is a vibrant orange with a pattern of concentric, overlapping circles.

Synergy Ventricular MENU DRIVEN ACCESS



Hummingbird™
Synergy Ventricular

MENU DRIVEN MONITORING

The Hummingbird Synergy Ventricular with FlexPort™ refines multi modal monitoring by providing a flexible access for your choice of probes. The menu driven approach provides concurrent monitoring of:

- ICP
- Ventricular Drainage
- FlexPort™ Probe Access



innerspace
Access for Life



Hummingbird™
Synergy Ventricular

FlexPort Probe Access

- Provides dedicated probe access with dimensions of 1.5mm diameter (.06") with a 1.2cm radius curve and a 30° exit angle
- Probe exits approximately 30mm below the dura
- Probe is directed into undisturbed tissue

AirPulse™ ICP

The Hummingbird Synergy brings operational **FLEXIBILITY** and **SIMPLICITY** to ICP monitoring utilizing the AirPulse™ Air Management System (AMS).

AirPulse AMS enhances ICP monitoring:

- Intra parenchymal sensing avoids loss of ICP
- Provides concurrent ventricular drainage and ICP monitoring
- **No** Capital Instrument
- **No** Leveling Required
- **Re-Zero** in Situ

AirPulse AMS

- Integrated on patient cable



MRI Conditional

The Hummingbird Synergy Ventricular is MRI Conditional, providing more flexible patient care and management.

innerspace

Access for Life

InnerSpace, Inc.

1622 Edinger Avenue, Suite C, Tustin, CA 92780
Customer Service: 877.486.2473
Customer Service Fax: 877-235-6902
www.innerspacemedical.com

Product	Model Number	Description	Kit Components	Order Quantity
Hummingbird Synergy Ventricular	HUMV-400MR	Monitoring kit for ICP and Ventricular Drainage. Multi-Parameter port for Oxygen or Blood Flow Probe	Drill Assembly Bolt Assembly Insert Assembly Catheter Assembly Bolt Stops Tripod PASS	4/Box

